

1. SET-UP

- a. Locate and turn on all EKG Patch Solution components; ensure the HD+ transmitter and touchECG application on your tablet are connected.
- b. Set-up the Patient Profile in the TouchECG application on your tablet.
- c. Place patient in a reclined position; ensure skin is clean and dry.

Don't forget to fill in all mandatory fields in the Patient Profile- especially Date of Birth!

2. ACQUIRE

- a. Select the correct size EKG Patch.
- b. Apply the EKG Patch and connect to the HD+ transmitter with the cable.
- c. Acquire the EKG using the touchECG application.

3. INTERPRET

- a. The EKG will automatically transmit for over-read services.
- b. Review the Unconfirmed EKG tracing to effectively triage the patient.
- c. Remove and dispose of the EKG Patch.
- d. Clean the HD+ transmitter and cable; return all components to the EKG Patch Solution Kit. Plug in the tablet to charge as needed.

